

ZEN (blue edition) 3.1


Guided Acquisition



Dr. Marion Lang
Product Management
2019-11-04

Guided Acquisition

Workflow Diagram




Guided Acquisition -Version 1.4

Guided Acquisition – Find objects of interest

Set-up experiments



1. Overview Experiment (e.g. tiles scan, low magnification)
Make sure to include the Objective in the Experiment setting!
Imaging Setup -> click on the objective -> activate „Objective Position“
2. Detailed Experiment (e.g. z-stack, higher magnification)
Make sure to include the Objective in the Experiment setting!




Guided Acquisition – Find objects of interest

Set-up image analysis



Execute the overview experiment manually (or a few tiles) to set up the image analysis:

- Detect the objects of interest
- Make sure not to detect too many objects for testing purposes
- Define measurement parameters for the objects (necessary for the script!):
 - Object ID
 - Bound Center Stage X
 - Bound Center Stage Y
 - Bound Width
 - Bound Height
 - ImageSceneContainerName



Guided Acquisition – Find objects of interest

Create a User Interface



```
# Initialize Dialog
window = ZenWindow()
window.Title = 'MyDialog'
window.Initialize('Guided Acquisition - Version : 1')

# add components to dialog
window.AddLabel('Image')
window.AddImage2dView(image)
window.AddDimension(image)
window.AddFolderBrowser('destfolder', 'Destination folder', 'C:\\\\OAD\\\\Output\\\\CZI Images')
window.AddDropDown('ddval', 'Dropdown', ['first', 'second', 'third'], 1)

# show the window
result = window.Show()
```

→ OAD_Training_User dialog with automatic arrangement

Guided Acquisition – Find objects of interest


Read out the Info from the User Interface



```
## get and check results
result=window.Show()

nosamples = str(result.GetValue('nosamples'))
concentration = str(result.GetValue('concentr'))

print(nosamples)
print(concentration)
```



Guided Acquisition – Find objects of interest

Generate a Drop-Down with Experiments and IAs




```
from System.IO import Directory, Path
```

```
def getshortfiles(filelist):  
    files_short = []  
    for short in filelist:  
        files_short.append(Path.GetFileName(short))  
  
    return files_short
```

```
# check the location of experiment setups and image analysis settings are stored  
docfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.UserDocuments)  
imgfolder = Zen.Application.Environment.GetFolderPath(ZenSpecialFolder.ImageAutoSave)
```

```
# get list with all existing experiments and image analysis setup and a short version of that list  
expfiles = Directory.GetFiles(Path.Combine(docfolder, 'Experiment Setups'), '*.czexp')  
ipfiles = Directory.GetFiles(Path.Combine(docfolder, 'Image Analysis Settings'), '*.czias')  
expfiles_short = getshortfiles(expfiles)  
ipfiles_short = getshortfiles(ipfiles)
```



→ OAD_Training_Load and show experiments and IA in dropdown

Guided Acquisition – Find objects of interest

Execute an Experiment



Set up an experiment, and test it in the GUI first.
Then execute via Script.

```
## Load experiment
exp.Load("OAD_Training_Exp", ZenSettingDirectory.User)
## Execute experiment, display acquired image
## Image is saved automatically in temp folder of AutoSavePath of Saving tab of Tools/Options menu
image = Zen.Acquisition.Execute(exp)
```

→ OAD_Training_Load and execute a tiles experiment


Guided Acquisition – Find objects of interest

Perform Autofocus - Simple



```
exp = ZenExperiment()  
exp = Zen.Acquisition.Experiments.ActiveExperiment  
#exp = Zen.Acquisition.Experiments.GetByName("testexp24.czexp")  
  
# get actual position  
old_pos = Zen.Devices.Focus.ActualPosition  
print(str(old_pos))  
  
# Run Autofocus based on the settings defined in the Experiment  
try:  
    Zen.Acquisition.FindAutofocus(exp, 30)  
except:  
    print("SWAF failed, the original position will be used")
```

→ OAD_Training_SWAF simple



Guided Acquisition – Find objects of interest

Perform Autofocus - Advanced



→ OAD_Training_SWAF_programmable

```
def runSWAF_special(SWAF_exp, delay=5, searchStrategy='Full',
                    sampling=ZenSoftwareAutofocusSampling.Coarse, relativeRangeIsAutomatic=False,
                    relativeRangeSize=500, timeout=0):

    # get current z-Position
    zSWAF = Zen.Devices.Focus.ActualPosition
    # set DetailScan active
    SWAF_exp.SetActive()
    # set SWAF parameters
    SWAF_exp.SetAutofocusParameters(searchStrategy=searchStrategy, sampling=sampling,
                                    relativeRangeIsAutomatic=relativeRangeIsAutomatic, relativeRangeSize=relativeRangeSize)
    SWAF_exp.Close()

    return zSWAF

exp = ZenExperiment()
exp = Zen.Acquisition.Experiments.GetByName("OAD_Training_Exp.czexp")
zSWAF_new = runSWAF_special(exp)
```

Guided Acquisition – Find objects of interest

Modify Tiles Region of an Experiment



```
exp = ZenExperiment()  
# Load a Tiles experiment!  
exp.Load("OAD_Training_exp.czexp")  
  
# get focus position  
zpos = Zen.Devices.Focus.ActualPosition  
  
# Perform Image Analysis and Create table with Bound and Stage information for each object  
SingleObj = Zen.Analyzing.CreateRegionTable(img)  
# get Bounds info from IA  
soi = SingleObj.GetBoundsColumnInfoFromImageAnalysis(True)  
  
# check the number of objects = rows inside image analysis table  
num_POI = SingleObj.RowCount
```

→ OAD_Training_Modify tile regions

Guided Acquisition – Find objects of interest

Modify Tiles Region of an Experiment



```
# execute experiment at the position of every detected object
```

```
for i in range(0, num_POI, 1):
```

```
    # get XY-stage position from center position of object
```

```
    xpos = SingleObj.GetValue(i, soi.CenterXColumnIndex)
```

```
    ypos = SingleObj.GetValue(i, soi.CenterYColumnIndex)
```

```
    # get height and width of bounding box of object
```

```
    bcwidth = SingleObj.GetValue(i, soi.WidthColumnIndex)
```

```
    bcheight = SingleObj.GetValue(i, soi.HeightColumnIndex)
```

```
    # Modify the XYZ position and size of the TileRegion on-the-fly
```

```
    exp.ClearTileRegionsAndPositions(0)
```

```
    exp.AddRectangleTileRegion(0, xpos, ypos, bcwidth, bcheight, zpos)
```

```
    # execute the tiles experiment
```

```
    output = Zen.Acquisition.Execute(exp)
```

Guided Acquisition – Find objects of interest

OAD Script



Get **Guided_Acquisition_shortUI.py** from https://github.com/zeiss-microscopy/OAD/tree/master/Guided_Acquisition

Try Guided Acquisition Macro and adapt (e.g. perform Image processing steps,...)

→ Guided_Acquisition_shortUI

